


# FIRST AID



SAFETY AT WORK

Eye Wash  
Plasters  
First Aid Stations

plum

# Optimize safety. Systematize first aid

With Plum's first aid concepts you will limit the extent and impact of accidents to the gratification of your employees and your company

---



Minor and major accidents happen every day. It is therefore important to have first aid at hand – right there where the accident happens. We offer you and your colleagues a wide range of first aid products so that you will be able to quickly treat cuts and minor wounds and eye injuries, safely and hygienically. We have also developed both first aid stations and mobile solutions that give an overview and ensure that first aid is always at hand.

## **CUSTOMIZED SOLUTION FOR QUICK AND EFFECTIVE HELP**

The first aid products include a wide range of eye wash solutions, plasters, wound cleansing wipes, etc. You can choose exactly what products you want. We have solutions for you whether you work in an office or handle acids/alkali in a high-risk environment, for example, when charging heavy batteries, industrial cleaning, in a laboratory or in the petrochemical industry.

## **” LOST MOST OF HIS EYESIGHT WHEN THE HOSE CAME LOOSE**

Hanstholm Fiskemelsfabrik has an advanced production facility, but a couple of years ago it only had ordinary eye wash at the plant. That cost one of the smiths most of his eyesight when a hose came loose and he got lye into both eyes. As a result of the accident, the factory upgraded its first aid solution:

“If we had had pH Neutral DUO then, the outcome of the accident would have been different. We have it now,” says Ronni Jørgensen, Head of Purchasing.



The green Plum Eye Wash helps you remove dirt and foreign objects

The blue pH Neutral provides immediate first aid in connection with serious accidents involving acids and alkali

The soft bottles are easy to handle and rinse the eye with a soft and even flow

There are no doubts about how to use the bottle. Clear pictograms show you how

The economic soft plastic eye cup is shaped to fit the curve of the eye. Squeeze it together before putting it to your eye. The cup will help keep the eye open and provide the best possible flow of liquid and rinsing time

Quick and easy to open. A single turn breaks the seal and you are ready to rinse

Economic and a long shelf-life of 3 years. The expiry date is clearly shown on the front of the label

With the ergonomic DUO eye cup, both eyes can be rinsed at the same time. The cup will help keep the eyes open and provide the best possible flow of liquid and rinsing time



## We cover the entire spectrum

The green Plum Eye Wash deals with dirt. The blue pH Neutral is crucial in connection with accidents involving acids and alkali

### YOU ONLY HAVE TWO EYES – TAKE CARE OF THEM

Your eyesight is one of the most important senses. So take care of your eyes. Regardless how careful you are accidents can happen, but everybody can make sure that they are prepared. With our eye wash solutions, you have done your best.

If an accident happens, seconds can decide whether your eyesight is damaged or not. Therefore, quick access to eye wash is crucial for all workplaces where there is a risk of eye injuries.

### A WIDE RANGE OF UNIQUE RINSING SOLUTIONS

Our wide product range ensures the best possible solution specifically for your needs. Choose small and large wall-mounted stations directly at the workplace. Or small, mobile solutions for the van, the tool box and the first aid box, and on the cleaning trolley in the special belt bag.



The bottles are CE marked according to the Guidelines for Classification of Medical Devices





# How to prevent eye injuries

Plum Eye Wash removes foreign objects like dust and dirt from your eye, while pH Neutral will neutralize acids and alkali

Our concepts have been developed in close collaboration with skilled eye specialists. This means that the products have been thought out down to the smallest detail. They are easy to use. Clear pictograms show you the correct use. The eye cup helps keep the eye open. The rinsing liquid hits the eye in the right place and with the right amount in a soft and even flow.



Plum Eye Wash is easy to use if an accident happens



Turn the eye cup in the direction of the arrow until the seal breaks



Lean back, squeeze the bottle and rinse



The eye can also be rinsed while you lean forward - by squeezing the bottle gently

## EYE WASH PLAN

1	2	3	4
Determine type of substance	Select appropriate eye wash and rinse	Continue rinsing with Plum Eye Wash	Rinse as directed
<b>Foreign substances</b> Dust, dirt, metal, wood, splinters, etc.	Plum Eye Wash	Plum Eye Wash	Rinse until substance is removed
<b>Acid</b>	pH Neutral	Plum Eye Wash	Rinse for about 2 minutes using pH Neutral and continue using Plum Eye Wash until you reach a doctor or an emergency room
<b>Alkali</b>	pH Neutral	Plum Eye Wash	Rinse for about 2 minutes using pH Neutral and continue using Plum Eye Wash until you reach a doctor or an emergency room
<b>Other chemicals</b>	Plum Eye Wash	Plum Eye Wash	Rinse until you reach a doctor or an emergency room

Always ensure the eye wash liquid hits the eye with a soft and even flow. In case of concurrent injury on both eyes, use DUO for rinsing both eyes simultaneously.



## Obvious and unique

The Plum DUO eye cup is a unique solution on the market. You no longer have to choose which eye to save

---

### **ABOUT 50% OF ALL ACCIDENTS AFFECT BOTH EYES**

Experience shows that about 50% of all eye accidents involving chemicals affect both eyes. With our DUO eye cup, you are prepared. You can rinse both eyes at the same time and will not have to choose which eye to save.

### **THE EYE CUP HELPS TO PROVIDE OPTIMAL RINSING**

Time is a critical factor, especially with accidents involving acids/alkali. Of course the injury depends on the strength of the acid or alkali, but apart from that, the contact time is crucial. The optimal first aid must therefore neutralize the harmful substance as quickly as possible. With the unique DUO cup, you will treat both eyes at the same time.

The eye cup is shaped to help you keep your eyes open during rinsing and to ensure a correct flow of liquid and rinsing time.

Plum Eye Wash and pH Neutral are disposable products in order to ensure that the rinsing liquid is sterile every time



## Eye wash in case of foreign objects

Plum Eye Wash can remove discomfort caused by foreign bodies in your eyes and prevent subsequent complications

---

It is easy to get a foreign body in your eye. Dust from the air or dirt, cleaning agents or paint from your hands. The risk is increased when you are working with machines and tools that can result in getting wood or metal splinters in your eyes.

Getting a foreign body in your eye can be extremely unpleasant, as the eye is extremely sensitive. Without treatment, it may develop into an inflammation of the eye, cornea or iris. By rinsing immediately after the incident, you will prevent this and prevent foreign bodies from getting stuck in the eye.

Plum Eye Wash contains a sterile saline solution of 0.9% sodium chloride, corresponding to the eye's natural fluid. The eye will be rinsed with a soft and pleasant flow. The product range has everything from the small handy bottle that you can always carry with you, to wall-mounted stations that will make it easy to get help in case of an emergency.



**200 ML  
PLUM EYE WASH**

A small and handy pocket-size bottle with 0.9% sterile sodium chloride. The bottle is easy to store in the specially designed belt bag, in the first aid kit, in the tool box, etc. It comes with a dust cap.

**Rinsing time:** about 2 mins.  
**Packaging:** 10 x 200 ml  
**Article no:** 4691



**500 ML  
PLUM EYE WASH**

A handy bottle with 0.9% sterile sodium chloride that can be used individually or as refill at eye wash stations and boxes. It comes with a dust cap.

**Rinsing time:** about 5 mins.  
**Packaging:** 12 x 500 ml  
**Article no:** 4604



**1000 ML  
PLUM EYE WASH**

Bottle with 0.9% sterile sodium chloride. Contains 1000 ml eye wash, which provides a longer rinsing effect. The bottle can be used individually or as refill at eye wash stations and boxes. It comes with a dust cap.

**Rinsing time:** about 10 mins.  
**Packaging:** 6 x 1000 ml  
**Article no:** 4707



**500 ML  
PLUM EYE WASH DUO**

A handy bottle with 0.9% sterile sodium chloride. The ergonomic DUO eye cup makes it possible to rinse both eyes at the same time. The bottle can be used individually or as refill at eye wash stations and boxes.

**Rinsing time:** about 2 mins.  
**Packaging:** 6 x 500 ml  
**Article no:** 4861



**1000 ML  
PLUM EYE WASH DUO**

Bottle with 0.9% sterile sodium chloride. Contains 1000 ml eye wash, which provides a longer rinsing effect. The ergonomic DUO eye cup makes it possible to rinse both eyes at the same time. The bottle can be used individually or as refill at DUO eye wash stations and boxes.

**Rinsing time:** about 5 mins.  
**Packaging:** 6 x 1000 ml  
**Article no:** 4800



**3 X 500 ML  
PLUM EYE WASH**

Practical and handy set with 3 x Plum Eye Wash bottles with 0.9% sterile sodium chloride. The box is easy to take with you in the car on the way to a doctor or the emergency room. The bottles can be used individually or as refills at eye wash stations and boxes.

**Rinsing time:** about 15 mins.  
**Packaging:** 4 x 3 x 500 ml  
**Article no:** 4620



# Eye wash in case of chemical accidents

pH Neutral neutralizes  
both acids and alkali.  
It also flushes them from  
the eye during rinsing.  
Proven safety and effect!

## IMPORTANT

Reservation is made for hydrofluoric acid, which is toxic and highly corrosive to skin and bones.

pH Neutral and Plum Eye Wash may neutralize and rinse, but not limit the corrosive effect, and therefore they are not recommended as first aid in connection with eye accidents involving hydrofluoric acid.

## RECOMMENDED USE OF PH NEUTRAL

pH Neutral is a two minutes eye wash solution for accidents involving acids and alkali. Rinse with pH Neutral in order to reduce the pH value of the harmful substance to a harmless level as quickly as possible. Then continue rinsing with Plum Eye Wash 0.9% sodium chloride until you reach a doctor or the emergency room.

## DANGEROUS CHEMICALS ARE EVERYWHERE

In a large number of jobs, businesses and industries, you are particularly exposed to serious eye injuries caused by chemicals. In airports, chemical factories, the recycling industry, in cement factories and on construction sites, in the food, metal, paper and mining industry, or in the petrochemical industry. Strong acids and alkali are used in the production and in connection with cleaning. The risk that something will go wrong, for example, because a hose comes loose or a battery explodes, is part of everyday life. You need pH Neutral at hand.

## JUST A FEW SECONDS CAN MAKE ALL THE DIFFERENCE

Acids and alkali break down the structure of skin and mucosa membranes. The severity of corrosion injuries depends primarily on:

- The concentration of the chemical.  
Strong acids and alkali are extremely corrosive.
- The temperature.  
A high temperature increases the corrosive properties of the chemical.
- The duration of contact.  
The longer the duration, the deeper the harmful substance will penetrate.

In connection with eye injuries involving acids and alkali, it is important to neutralize the chemical as quickly as possible to a pH level of about 7.4, corresponding to the normal pH level of the eye's fluid - preferably before the acid or alkali has penetrated into the cornea. To optimize safety where there is a risk of acid and alkali injuries, the eye wash solution must be:


- Easily accessible
- Quick to use
- Effective

## THE EYE WASH LIQUID MUST BE EFFECTIVE

pH Neutral is a 4.9% phosphate buffer that will quickly and effectively neutralize even very strong acids and alkali. The liquid consists of phosphate salts that are present in the body. The hydrogen phosphate buffer is effective against acids, whereas the dihydrogen phosphate buffer reacts against alkali. The effect of pH Neutral has been documented in scientific studies.<sup>1) 2)</sup>

We have tested the neutralising effect of pH Neutral for both strong acids and alkali compared with the effect of ordinary eye wash with 0.9% sodium chloride. The test clearly shows that pH Neutral will quickly neutralize both acids and alkali to a harmless pH level. Ordinary eye wash dilutes the liquid, but only small changes to the pH value are obtained (Figure 1).

Figure 1: Neutralization of strong alkali and strong acid



1: Prof. Dr. med. Schrage N.F., Facharzt für Augenheilkunde an der Augenklinik Merheim, Ostmerheimer Str. 200 in 51109 Köln, in Kooperation mit dem Aachener Centrum für Technologietransfer in der Ophthalmologie ACTO, Geschäftsstelle Karlsburgweg 9 in 52070 Aachen, Vortrag: Augenverätzungen. Was gibt es Neues? Regionalforum Arbeitsmedizin, Stuttgart 23.02.2006, Folien 44 ff.

2: Prof. Dr. med. Schrage N.F., Facharzt für Augenheilkunde an der Augenklinik Merheim, Ostmerheimer Str. 200 in 51109 Köln, in Kooperation mit dem Aachener Centrum für Technologietransfer in der Ophthalmologie ACTO, Geschäftsstelle Karlsburgweg 9 in 52070 Aachen, Vergleichsergebnisse von pH Neutral an verätzten Augen, Spülversuche an mit NaOH verätzten Kanichenaugen im Vergleich von Plum pH Neutral mit physiologischer Kochsalzlösung, 24.05.2005.

All documentation is available on request.





**200 ML  
PH NEUTRAL**

A small, handy pocket-size bottle with 4.9% sterile phosphate buffer that will quickly neutralize acids and alkali. The bottle is easy to carry in the specially designed belt bag, in the first aid box, in the tool box, etc. It can also be used as refill at eye wash stations and boxes. It comes with a dust cap.

**Rinsing time:** about 2 mins.  
**Packaging:** 10 x 200 ml  
**Article no:** 4752



**500 ML  
PH NEUTRAL DUO**

A handy bottle with 4.9% sterile phosphate buffer that will quickly neutralize acids and alkali. The ergonomic DUO eye cup makes it possible to rinse both eyes at the same time. It is particularly suitable for workplaces where there is a risk that both eyes are affected. The bottle can be used individually or as refill at eye wash stations and boxes.

**Rinsing time:** about 2 mins.  
**Packaging:** 6 x 500 ml  
**Article no:** 4801



**1000 ML  
PH NEUTRAL EYEWASH SHOWER**

A bottle with 4.9% sterile phosphate buffer that will quickly neutralize acids and alkali. The bottle can be used individually at the workbench, in the tool box, the van, or mounted on the wall close to the workplace. It comes with a dust cap.

**Rinsing time:** about 2 mins.  
**Packaging:** 6 x 1000 ml  
**Article no:** 4746



**” ACCIDENTS INVOLVING  
ACIDS AND ALKALI WILL  
HAPPEN WHEN YOU  
LEAST EXPECT THEM**

pH Neutral DUO saved the eyesight of a SITA France employee. They handle plastic waste such as plastic bottles with cleaning detergents. The liquid from a cleaning bottle went into both the employee's eyes. Thanks to pH Neutral DUO, the harmful substance was neutralized, and the employee was already back at work the next day.





# Eye wash stations for accidents involving foreign bodies

A wall-mounted solution gives  
a good overview and help is right at hand

## SHOW THAT YOU TAKE FIRST AID SERIOUSLY

When an accident has happened, it is too late to look for help in a box or under a cupboard.

With the wall-mounted stations, you will systematize first aid and make it easy to find when the accident happens. Choose between open wall stations and box solutions. The latter are particularly suitable for dusty work environments in order to protect the bottles and avoid adding more dust to an accident.

Always have additional bottles at hand so that rinsing may continue on the way to further treatment.

### PLUM EYE WASH BOX DUO

Dust-proof wall station in impact-proof polystyrene with 2 x 1000 ml Plum Eye Wash. Particularly suitable for small and mobile workplaces with a large risk of getting dust and dirt into both eyes.

**Rinsing time:** about 10 mins.  
**Packaging:** 1  
**Article no:** 4816  
**Dimensions:** h: 357 x w: 252 x d: 116 mm



### PLUM EYE WASH SINGLE CARTON

Single carton with 200 ml Plum Eye Wash. The carton protects the bottle and is also easy to open when it is needed. Suitable for refilling stations and for mobile workplaces, such as fitters' vans and goods transport.

**Rinsing time:** about 2 mins.  
**Packaging:** 10  
**Article no:** 4695  
**Wall bracket art. no:** 955904  
**Dimensions:** h: 58 x w: 58 x d: 163 mm



### PLUM EYE WASH STATION SINGLE CARTON

Single carton with 500 ml Plum Eye Wash and wall bracket. The carton protects the bottle and is also easy to open when it is needed. Suitable for quick and easy eye rinsing in small workplaces with dust and dirt.

**Rinsing time:** about 5 mins.  
**Packaging:** 12  
**Article no:** 4609  
**Refill article no:** 4605  
**Wall bracket art. no:** 955904  
**Dimensions:** h: 79 x w: 79 x d: 214 mm



### PLUM EYE WASH STATION DUO

Wall station with 1000 ml Plum Eye Wash DUO. Particularly suitable for small and mobile workplaces where prolonged rinsing of both eyes may be required.

**Rinsing time:** about 5 mins.  
**Packaging:** 1  
**Article no:** 4802  
**Dimensions:** h: 375 x w: 180 x d: 95 mm



### PLUM EYE WASH BOX

Dust-proof wall station in impact-proof polystyrene with 2 x 500 ml Plum Eye Wash. Particularly suitable for production areas with a lot of dust and dirt, and for mobile workplaces.

**Rinsing time:** about 10 mins.  
**Packaging:** 1  
**Article no:** 4650  
**Dimensions:** h: 280 x w: 230 x d: 110 mm

### ALSO AVAILABLE WITH A HEATING UNIT

**Article no:** 4649

### EASY TO MOUNT

All wall stations are easy to mount and come with pictograms that clearly indicate how to use the products.



### PLUM EYE WASH STATION

Wall station with 500 ml Plum Eye Wash. Particularly suitable for small and mobile workplaces.

**Rinsing time:** about 5 mins.  
**Packaging:** 1  
**Article no:** 4611  
**Dimensions:** h: 300 x w: 170 x d: 80 mm



### PLUM EYE WASH STATION MAXI

Wall station with 2 x 1000 ml Plum Eye Wash. Particularly suitable for places where prolonged rinsing may be required. It comes with a dust cap.

**Rinsing time:** about 20 mins.  
**Packaging:** 1  
**Article no:** 4708  
**Dimensions:** h: 375 x w: 260 x d: 95 mm



### PLUM EYE WASH STATION

Wall station with 2 x 500 ml Plum Eye Wash. Particularly suitable for environments and small workplaces that are not dirty, but where there is a risk of getting foreign bodies into the eyes.

**Rinsing time:** about 10 mins.  
**Packaging:** 1  
**Article no:** 4694  
**Dimensions:** h: 290 x w: 228 x d: 80 mm



# Eye wash stations for accidents with acids and alkali

With a combi-box, you are covered  
- from minor incidents to serious accidents involving acids and alkali

Suddenly a pipe with liquid under pressure bursts - or a hose comes off a connector. As a result, heavily corrosive liquids spurt out. That is how the most frequent accidents with chemical liquids happen. If you get it into your eyes, every second counts if your eyesight is to be saved.

## PLUM EYE WASH AND PH NEUTRAL IS A STRONG TEAM

In case of accidents involving acids and alkali, we recommend a combined rinsing – first with pH Neutral, then with Plum Eye Wash. We have therefore made a selection of combi-boxes so that you can start by neutralising the acid or alkali, and then continue to rinse until you reach a doctor or the emergency room.

If you work mostly with liquid chemicals, we recommend pH Neutral DUO, so that you can rinse both eyes in case of an emergency.

We recommend the open stations for environments without dust and dirt, such as for example laboratories. The closed stations protect the products in dusty environments. All stations are easy to mount and come with a mirror and pictograms clearly indicating how to use the eye wash bottles.



### PH NEUTRAL SINGLE CARTON

Single carton with 200 ml pH Neutral. The robust carton protects the bottle and is also easy to open when it is needed. For quick neutralization of acids and alkali. Suitable for refilling stations and for mobile workplaces, such as fitters' vans and goods transport.

**Rinsing time:** about 2 mins.  
**Packaging:** 10  
**Article no:** 4755  
**Wall bracket art. no:** 955904  
**Dimensions:** h: 58 x w: 58 x d: 163 mm



### QUICKCLEAN WOUND CLEANSING WITH WALL BRACKET

Box with 40 wound cleansing wipes for cleaning minor cuts and wounds. Comes with a wall bracket and can be placed as close to the workplace as possible. Each wipe is packed separately and contains sterile water. The wipes measure 130 x 200 mm.

**Packaging:** 6 x 40  
**Article no:** 5550  
**Dimensions:** h: 150 x w: 80 x d: 64 mm



### FIRST AID STATION

Wall station with Plum Eye Wash, QuickFix plaster dispenser, QuickClean wound cleansing wipes, mirror and pictogram.

**CONTENTS:**  
1 x 500 ml Plum Eye Wash 4604  
1 x 40 QuickClean wound cleansing wipes 5550  
1 x QuickFix plaster dispenser with 90 plasters (Elastic & Water Resistant) 5507

**Packaging:** 1  
**Article no:** 5101  
**Dimensions:** h: 300 x w: 515 x d: 85 mm





### PLUM COMBI-BOX DUO

Dust-proof wall station in impact-proof polystyrene with 500 ml pH Neutral DUO and 500 ml Plum Eye Wash DUO. Particularly suitable for dusty production areas where there is a risk of accidents involving foreign bodies and acids and alkali, and where there is a risk that both eyes will be affected.

**Rinsing time:** pH Neutral about 2 mins.  
Plum Eye Wash about 2 mins.  
**Packaging:** 1  
**Article no:** 4862  
**Dimensions:** h: 280 x w: 230 x d: 110 mm



### PLUM COMBI-BOX DUO MAXI

Dust-proof wall station in impact-proof polystyrene with 500 ml pH Neutral DUO and 1000 ml Plum Eye Wash DUO. Particularly suitable for dusty production areas where there is a risk of accidents involving foreign bodies and acids and alkali, and where there is a risk that both eyes will be affected.

**Rinsing time:** pH Neutral about 2 mins.  
Plum Eye Wash about 5 mins.  
**Packaging:** 1  
**Article no:** 4810  
**Dimensions:** h: 357 x w: 252 x d: 116 mm



### PLUM COMBI-BOX

Dust-proof wall station in impact-proof polystyrene with 200 ml pH Neutral and 500 ml Plum Eye Wash. Particularly suitable for dusty production areas where there is a risk of accidents involving foreign bodies and acids and alkali.

**Rinsing time:** pH Neutral about 2 mins.  
Plum Eye Wash about 5 mins.  
**Packaging:** 1  
**Article no:** 4787  
**Dimensions:** h: 280 x w: 230 x d: 110 mm

**ALSO AVAILABLE WITH A HEATING UNIT**  
**Article no:** 4664



### PH NEUTRAL EYEWASH SHOWER STATION

Wall station with 1000 ml pH Neutral, which will quickly neutralize acids and alkali. Suitable for mounting on the wall near the workplace. It comes with a dust cap.

**Rinsing time:** about 2 mins.  
**Packaging:** 1  
**Article no:** 4741  
**Dimensions:** h: 300 x w: 160 x d: 95 mm



### PLUM COMBI-STATION

Wall station with 200 ml pH Neutral and 500 ml Plum Eye Wash. Particularly suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eye.

**Rinsing time:** pH Neutral about 2 mins.  
Plum Eye Wash about 5 mins.  
**Packaging:** 1  
**Article no:** 4773  
**Dimensions:** h: 290 x w: 228 x d: 80 mm



### PLUM COMBI-STATION DUO

Wall station with 500 ml pH Neutral DUO and 1000 ml Plum Eye Wash DUO. Particularly suitable for places where there is a risk of getting both foreign bodies and acids and alkali in the eyes and where there is a risk that both eyes will be affected.

**Rinsing time:** pH Neutral about 2 mins.  
Plum Eye Wash about 5 mins.  
**Packaging:** 1  
**Article no:** 4803  
**Dimensions:** h: 375 x w: 250 x d: 95 mm

# Quick

# FIX

## QUICKFIX IS A UNIQUE DISPENSER:

- Gives immediate and visible access to plasters and eye wash ampoules
- Each plaster is ready for use just by pulling it from the dispenser
- Each plaster is packed separately
- Only one hand is needed to apply
- Can be refilled as required
- Easy to see when it is empty
- Easy to stop minor bleeding, also without staining the other plasters
- Small dimensions - easy to find room for the QuickFix dispenser
- Includes a shelf for storing small, relevant objects
- Can be locked so that it can only be opened with a key - may prevent unwanted waste



### QUICKFIX ELASTIC

Wall-mounted plaster dispenser with 90 elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.

<b>Packaging:</b>	1
<b>Article no:</b>	5502
<b>Refill article no:</b>	5512
<b>Dispenser dimensions:</b>	h: 135 x w: 230 x d: 32 mm
<b>Plaster dimensions:</b>	l: 72 x w: 25 mm



### QUICKFIX WATER RESISTANT

Wall-mounted plaster dispenser with 90 water resistant plasters that are particularly suitable for working environments with a high level of humidity. The plasters are made of a protective PE material, which allows the skin to breathe.

<b>Packaging:</b>	1
<b>Article no:</b>	5501
<b>Refill article no:</b>	5511
<b>Dispenser dimensions:</b>	h: 135 x w: 230 x d: 32 mm
<b>Plaster dimensions:</b>	l: 72 x w: 25 mm



### QUICKFIX ELASTIC & WATER RESISTANT

Wall-mounted plaster dispenser with 45 elastic and 45 water resistant plasters. A combination of plasters covering many requirements at the workplace.

<b>Packaging:</b>	1
<b>Article no:</b>	5507
<b>Refill article no:</b>	5512 & 5511
<b>Dispenser dimensions:</b>	h: 135 x w: 230 x d: 32 mm
<b>Plaster dimensions:</b>	l: 72 x w: 25 mm

# QuickFix is quick help for minor day-to-day accidents

Mount the plaster on the wall and save time.  
Put the plaster on and back to work!

## QUICK AND HYGIENIC HELP WITH MINOR INJURIES

Lots of minor accidents happen every day in the workplace. We are talking about minor events where a plaster is all that is needed. Choose the right plaster for your work situation. Or the right combination. We have plasters for ordinary work situations, for workplaces with a lot of humidity and for work within the food sector. Do not waste time on minor injuries. Take a plaster from the dispenser, put it on, and go back to work! QuickFix solves the problem easily, hygienically and effectively.



Take the plaster from the dispenser.  
It is ready for use...



... and can be put on  
using one hand



### QUICKFIX DETECTABLE

Wall-mounted plaster dispenser with 90 blue and conspicuous plasters with built-in metal surface that can be detected by a metal detector. Particularly suitable for the food industry. The plasters are made from a breathable material.

**Packaging:** 1  
**Article no:** 5503  
**Refill article no:** 5513  
**Dispenser**  
**dimensions:** h: 135 x w: 230 x d: 32 mm  
**Plaster dimensions:** l: 72 x w: 25 mm



### QUICKFIX ELASTIC LONG

Wall-mounted plaster dispenser with 60 long and elastic textile plasters that adjust to the movement of the skin. The plasters are skin-friendly and allow the skin to breathe. The long plaster can be wrapped around the finger several times and is more firmly attached.

**Packaging:** 1  
**Article no:** 5528  
**Refill article no:** 5508  
**Dispenser**  
**dimensions:** h: 135 x w: 230 x d: 32 mm  
**Plaster dimensions:** l: 120 x w: 20 mm

## CHOOSE BETWEEN THREE TYPES OF PLASTER AND TWO LENGTHS

### QUICKFIX ELASTIC

Skin-friendly and comfortable plaster. Elastic so that it adjusts to the movements of the skin and the hand.

### QUICKFIX WATER RESISTANT

Waterproof plaster for working environments with a high level of humidity.

### QUICKFIX DETECTABLE

Special plaster for the food industry. The blue color and the built-in metal surface make it easy to trace the plaster both visually and with a metal detector.

### QUICKFIX LONG

Extra-long plasters. They go round the finger several times, which means that they are firmly attached.



### QUICKFIX DETECTABLE LONG

Wall-mounted plaster dispenser with 60 long blue and conspicuous plasters with built-in metal surface which can be traced by a metal detector. Particularly suitable for the food industry. The long plaster can be wrapped around the finger several times and is more firmly attached.

**Packaging:** 1  
**Article no:** 5529  
**Refill article no:** 5509  
**Dispenser**  
**dimensions:** h: 135 x w: 230 x d: 32 mm  
**Plaster dimensions:** l: 120 x w: 20 mm



#### QUICKFIX&RINSE

Wall-mounted dispenser with 5 eye wash ampoules with 0.9% sterile sodium chloride and 45 breathable and elastic textile plasters. The dispenser gives easy, quick and hygienic help in case of minor accidents with a possibility of combining plasters and eye wash as required.

---

<b>Packaging:</b>	1
<b>Article no:</b>	4631
<b>Refill article no:</b>	5512 & 5160
<b>Dispenser dimensions:</b>	h: 135 x w: 230 x d: 32 mm
<b>Plaster dimensions:</b>	l: 72 x w: 25 mm



## Small ampoules for small accidents

QuickRinse is an economical solution  
for small accidents



#### QUICKRINSE

Wall-mounted dispenser with 10 eye wash ampoules with 20 ml 0.9% sterile sodium chloride. The dispenser provides easy, quick and hygienic help, rinsing dust and dirt from the eyes. Or to regain comfort when you have dry or irritated eyes. Easy to mount where help is most needed.

---

<b>Packaging:</b>	1
<b>Article no:</b>	4632
<b>Refill article no:</b>	5160
<b>Dispenser dimensions:</b>	h: 135 x w: 230 x d: 32 mm

#### USE QUICKRINSE FOR MINOR RINSING

The eye is extremely sensitive. This means that dry eyes and minor irritation may lead to much discomfort with stinging, increased light-sensitivity and pain. Use QuickRinse and regain your well-being. Each ampoule contains 20 ml which will solve the problem caused by dry or irritated eyes. QuickFix&Rinse provides an optimum combination, so that you will have both plasters and eye wash ampoules at hand.





# Always at hand

## Practical solutions

– so that you can take first aid with you

Our plasters and small eye wash ampoules are also available in practical packaging that makes it easy for you to take the first aid with you. Place them in the glove compartment, the tool box, your bag, in a pocket, or an ordinary first aid box. Then you will be prepared for the minor accidents that always happen inconveniently and when you least expect it.



### QUICKFIX MINI

Small and practical pack with 30 elastic plasters. Easy to carry anywhere that you may need a plaster. Each plaster is packed separately and is ready for use just by pulling it from the pack. The plasters are skin-friendly and comfortable, as they adjust to your movements and allow the skin to breathe.

**Packaging:** 30 x 30  
**Article no:** 5504  
**Plaster dimensions:** l: 72 x w: 25 mm



### QUICKRINSE MINI

Small and practical pack with 4 eye wash ampoules, each with 20 ml 0.9% sterile sodium chloride. A handy solution that gives quick relief from dust and dirt in the eye or for eye-irritation caused by a dry environment. The ampoules open quick and easily with a single twist and are ready to use. The small pack size makes it easy to carry QuickRinse Mini with you.

**Packaging:** 5 x 4  
**Article no:** 4619



### BELT BAG FOR 200 ML BOTTLES

Practical bag for your belt, work clothes, etc. Easy to open and provides quick first aid in case of accidents. A safe and hygienic way to store and protect the bottle.

**Packaging:** 10  
**Article no:** 4692



# QuickSafe® gives you quick, easy and safe first aid

Get an overview and peace of mind that you will always have what you need in the right place

---

## CUSTOMIZED FIRST AID STATIONS


Fill your first aid station according to your needs. Eye wash for rinsing foreign bodies from the eye. pH Neutral for accidents with acids or alkali. Plasters that are suitable for your work situation. Burn gel that soothes minor burns. Eye wash ampoules for minor accidents and dry, irritated eyes. Wipes that will cleanse the wound. Compress dressings for more serious accidents. All of it behind a transparent cover that will give you a perfect overview of what you need.

## ” GET THE MOST IMPORTANT ITEMS IN THE RIGHT PLACE

Spæncom A/S is a leading supplier of concrete elements. The company has chosen Plum Quicksafe® in order to increase safety for their more than 500 employees. “We have chosen Plum Quicksafe®, because the station can be organized to suit the individual work area. In the areas where concrete elements are acid-washed, we can, for example, place Plum Eye Wash and pH Neutral DUO close to the work situation. QuickSafe® contains the most necessary items, it can be placed so that first aid is close at hand, and it is very practical with the transparent cover.”

Jany Johannsen,  
Environment and Safety Coordinator





### QuickSafe® EMPTY

**YOU MAKE THE DECISIONS**  
 If you want full freedom to organize your first aid according to your requirements, you can start with an empty station and fill it to suit your specific situation.



### BUY AN EMPTY STATION AND FILL IT TO MEET YOUR NEEDS

Of course you can get a 100% customized solution if you have requirements that we do not meet with the 4 different QuickSafe® solutions, which we present on the following page. You can start with a wall station and fill it with eye wash bottles and other first aid products.

**Packaging:** 1  
**Article no:** 5173  
**Dimensions:** h: 430 x w: 253 x d: 92 mm

# Help is at hand

Easy access to first aid is the key to increased safety

With QuickSafe® you gather the most important first aid products in one place and position the solution close to where the products are to be used. The station does not take up much more room than an A3 sheet of paper, so you can place QuickSafe® exactly where it is needed.

The transparent cover gives you an overview and provides quick and easy access to first aid. You can easily see when the station needs a refill, and the dust-proof cover together with the separately packed products make QuickSafe® both a hygienic and an elegant solution.



**SCAN THE CODE OR VISIT  
 PLUM-QUICKSAFE.COM  
 AND SEE HOW EASY IT IS TO  
 GET A FIRST AID SOLUTION**







## QuickSafe® BASIC

### A GOOD STARTING POINT

With a QuickSafe® Basic you can choose the contents for the top shelves. It gives you full flexibility.



### DO YOU HAVE SPECIAL REQUIREMENTS REGARDING YOUR FIRST AID STATION?

The backbone of the entire QuickSafe® series is Plum's acclaimed eye wash. With QuickSafe® Basic you will get it in a combination with pH Neutral, eye wash ampoules and plasters. The rest is up to you. Organize the first shelf in the station to suit your company's requirements.

CONTENTS:	Article no:
2 x 500 ml Plum Eye Wash	4604
1 x 200 ml pH Neutral	4752
5 x 20 ml QuickRinse eye wash ampoules	5160
1 x QuickFix Elastic refill with 45 plasters	5512

---

<b>Packaging:</b>	1
<b>Article no:</b>	5170
<b>Dimensions:</b>	h: 430 x w: 253 x d: 92 mm



## QuickSafe® COMPLETE

### THE COMPLETE SOLUTION

Eye wash including pH Neutral and eye wash ampoules, pressure dressing, wound cleansing wipes, plasters and burn gel to soothe minor burns.



### WIDE FIRST AID COVER

With QuickSafe® Complete you will get a wide range of first aid products. There are products for most critical situations and for many of the day-to-day minor accidents.

CONTENTS:	Article no:
18 x QuickCool burn gel	5150
1 x QuickStop wound compress dressing 1 x 17 x 17 cm and 2 x 12 x 7.5 cm	5152
20 x QuickClean wound cleanser	5151
2 x 500 ml Plum Eye Wash	4604
1 x 200 ml pH Neutral	4752
5 x 20 ml QuickRinse eye wash ampoules	5160
1 x QuickFix Elastic refill with 45 plasters	5512

---

<b>Packaging:</b>	1
<b>Article no:</b>	5174
<b>Dimensions:</b>	h: 430 x w: 253 x d: 92 mm





## QuickSafe® FOOD INDUSTRY

### FOR THE FOOD INDUSTRY

We have filled a station with the relevant content to suit the requirements of the food industry.



### QUICK AND SAFE HELP WITH BUILT-IN TRACEABILITY

A first aid station targeting the food industry, kitchens, canteens and restaurants with, among other things, traceable plasters, wound cleansing, burn gel and wound compress dressings.

CONTENTS:	Article no:
18 x QuickCool burn gel	5150
1 x QuickStop wound pressure dressing 1 x 17 x 17 cm and 2 x 12 x 7.5 cm	5152
20 x QuickClean wound cleanser	5151
2 x 500 ml Plum Eye Wash	4604
1 x 200 ml pH Neutral	4752
1 x QuickFix Detectable refill with 45 plasters	5513
1 x QuickFix Detectable Long refill with 30 plasters	5509

<b>Packaging:</b>	1
<b>Article no:</b>	5175
<b>Dimensions:</b>	h: 430 x w: 253 x d: 92 mm



## QuickSafe® CHEMICAL INDUSTRY

### STRONG SOLUTION FOR TOUGH PRODUCTS

Accidents with acids/alkali can have serious and life-long consequences for the individual if you compromise on accessibility.



### WHEN ONLY SECONDS STAND BETWEEN A MINOR AND A MAJOR ACCIDENT

If there is a risk of accidents involving acids or alkali, safety in the form of effective and mobile eye wash products is crucial. It contains, among other things, pH Neutral DUO that will rinse both eyes at the same time.

CONTENTS:	Article no:
1 x 500 ml Plum Eye Wash DUO	4861
1 x 500 ml pH Neutral DUO	4801
1 x 200 ml pH Neutral	4752
5 x 20 ml QuickRinse eye wash ampoules	5160
1 x QuickFix Elastic refill with 45 plasters	5512

<b>Packaging:</b>	1
<b>Article no:</b>	5171
<b>Dimensions:</b>	h: 430 x w: 253 x d: 92 mm



**QUICKSTOP**

Wound compress dressing set for effective first aid in the event of minor as well as major bleeding. Sterile dressing with integrated compression. Elastic, non-adhesive and skin-friendly. QuickStop for general use as pressure, support and protection dressing.

**THE SET INCLUDES:**

- 1 x wound compress dressing, 17 x 17 cm
- 2 x wound compress dressings, 12 x 7.5 cm

**Packaging:** 6 x 3  
**Article no:** 5152



**QUICKCLEAN WOUND CLEANSING**

Box with 40 wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Can be used as refill for the first aid station and QuickClean Wound Cleansing with wall bracket (p. 12). The wipes measure 130 x 200 mm

**Packaging:** 13 x 20  
**Article no:** 5551



**QUICKCLEAN**

Wound cleansing wipes for cleaning minor wounds and cuts. Each wipe is packed separately and contains sterile water. Can be used individually and for refilling the QuickSafe® stations. 20 wipes per box.

**Packaging:** 6 x 20  
**Article no:** 5151



**QUICKCOOL**

Burn gel for first aid in case of minor burns and scalding. Cooling and pain-relieving effect. Contains tea tree oil. Can be used individually and for refilling the QuickSafe® stations. 18 wipes per box.

**Packaging:** 6 x 18  
**Article no:** 5150

# List of refills

## Keep your first aid updated

You can easily keep your first aid stations and QuickFix dispensers filled and up-to-date. At the same time, you will be able to add new products to your first aid if the need arises. All stations have built-in flexibility so that you can put together a solution that suits your company.



### QUICKFIX WATER RESISTANT REFILL

45 water resistant plasters that are particularly suitable for working environments with a high level of humidity. The plasters are made from a protective PE material, which allows the skin to breathe.

**Packaging:** 6 x 45  
**Article no:** 5511  
**Plaster:** l: 72 x w: 25 mm



### QUICKFIX ELASTIC REFILL

45 elastic and comfortable textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe.

**Packaging:** 6 x 45  
**Article no:** 5512  
**Plaster:** l: 72 x w: 25 mm



### QUICKFIX DETECTABLE REFILL

45 blue and conspicuous plasters with a built-in metal surface that can be traced by a metal detector. Particularly suitable for the food industry. The plasters are made from a breathable material.

**Packaging:** 6 x 45  
**Article no:** 5513  
**Plaster:** l: 72 x w: 25 mm



### QUICKFIX ELASTIC LONG REFILL

30 long and elastic textile plasters that adjust to the movements of the skin. The plasters are skin-friendly and allow the skin to breathe. The long plasters can be wrapped around the finger several times and are more firmly attached.

**Packaging:** 6 x 30  
**Article no:** 5508  
**Plaster:** l: 120 x w: 20 mm



### QUICKFIX DETECTABLE LONG REFILL

30 long, blue and conspicuous plasters with built-in metal surface that can be traced by a metal detector. Particularly suitable for the food industry. The long plasters can be wrapped around the finger several times and are more firmly attached.

**Packaging:** 6 x 30  
**Article no:** 5509  
**Plaster:** l: 120 x w: 20 mm



### QUICKRINSE REFILL

5 eye wash ampoules with 20 ml 0.9% sterile sodium chloride. For refilling QuickSafe®, QuickFix&Rinse and QuickRinse dispensers.

**Packaging:** 4 x 5  
**Article no:** 5160



Plum A/S  
Frederik Plums Vej 2 · DK-5610 Assens  
Telefon +45 6471 2112  
info@plum.dk · www.plum.eu

plum